# **SonoTip<sup>®</sup> TopGain**

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**Warning !**

**For single use only. Do not reuse, reprocess or resterilize.** Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

**Caution:** Federal law restricts this device to sale by or on the order of a physician.

**Intended Purpose**

The SonoTip® TopGain needle system is intended for endoscopic ultrasound-guided fine needle biopsy (FNB) of the submucosal and extra-luminal lesions of the gastrointestinal tract (i.e., pancreatic masses, mediastinal masses, peri-pancreatic masses and lymph nodes).

**Contraindications**

The contraindications are those to the primary endoscopic procedure for providing access to the desired site.

The relative contraindications include, but not limited to: Coagulopathy.

**Possible Complications**

Possible complications from the extra-luminal endoscopic ultrasound-guided fine needle biopsy (FNB) may include but are not limited to infection, hemorrhage, perforation and tumor seeding. Extra-luminal FNB of cystic lesions have a higher risk of complication from infection!

**Warning!**

The user of this equipment should be thoroughly trained in the applicable procedure. Furthermore, failure to read and thoroughly understand the contents of this instruction manual may result in serious injury to the patient and/or user. It is essential to follow the instructions enclosed, as well as instructions found in other instruction manuals, infection control requirements and current hospital protocols for usage, cleaning and sterilization. Possible injuries related to endoscopic procedures may include perforation, electrical burns and shock, hemorrhage, infection, explosion, etc.

! Failure to follow these instructions may result in malfunction and/or damage of the instrument or endoscope!

**Storage**

Please store the instrument at room temperature, keep dry and avoid exposure to sunlight.

**General Safety Instruction**

- Wearing protective garments (i.e., gloves, protective eyewear, gowns, etc.) will assist in the protection against risk of cross-contamination and injury.
- Do not use an instrument if it does not function properly. If the device has any physical defect(s), replace the instrument with a new and/or functioning instrument.
- Check the endoscope's working channel diameter prior to introducing the instrument into the endoscope. To ensure the appropriate fit of the instrument the diameter of the endoscope's working channel needs to be at least "0.2 mm" larger than the instrument's outer diameter.

- ! Never use products which do not comply with the technical specifications!
- Do not force an instrument down the working channel of the endoscope. This could result in endoscope and/or instrument damage.
- Do not under any circumstances make any modifications to the instrument or components described in this manual.

### **Inspection and Testing**


- Note any damage/open seal of the packaging.
- If damage or an open seal is detected, sterilization may be compromised and the instrument should be discarded and replaced.
- Check date of sterilization.
- Never use products with non-valid sterilization date.
- Check the instrument prior to each use.
- ! Should any irregularities occur, dispose the instrument and replace it by a new one.
- ! Using a damaged instrument with physical or functional abnormalities may result in endoscope damage, patient and/or user injury.

### **Preparation of the Aspiration Syringe**


- Open package and remove the sterile aspiration syringe.
- Check to ensure that the stopcock valve is securely attached to the syringe.
- Ensure syringe plunger is fully encased in the syringe chamber.
- Turn the stopcock valve to its closed position.
- Pull the plunger completely back. Turn the plunger clockwise until the plunger is locked, thereby creating negative suction pressure.
- Set the prepared syringe aside until you are ready to use it for biopsy.

### **Preparation of the SonoTip® TopGain Needle**

- Remove the SonoTip® TopGain needle from its packaging and protective sleeve by unlocking from the Luer-Lock at the distal end of the handle (9) and sliding off the plastic sleeve.
- Inspect the Plastic Sheath (10) and verify that it is free of any noticeable kinks and/or defects and is fixed properly to the handle (8).
- ! Verify that the Sliding Sheath Handle (8) is fully retracted. **Picture 1** shows the correct position on the Scale for proper Sheath Adjustment (6).

- ! If it is determined that the Sliding Sheath Handle (8) is not fixed in the correct position as shown in **Picture 1**, please unlock the blue Sheath Adjuster (7) by turning the Sheath Adjuster (7) counter-clockwise. Retract the Sliding Sheath Handle (8) until it stops as illustrated in **Picture 1**. Turn the Sheath Adjuster (7) clockwise to lock the Sliding Sheath Handle (8) in position.
- ! Ensure that the Needle (11) and the Nitinol Stylet (1) are securely retracted into the Plastic Sheath (10) before advancing the instrument into the accessory channel of the endoscope. **Picture 2** shows the correct position on the Scale Needle Adjustment (5).

**Picture 2**

- ! If it is determined that the Scale for Adjustment (5) is not in the correct position as shown in **Picture 2**, retract the Needle Handle (2) fully until it stops and is in position as illustrated in **Picture 2**. Unlock the Needle Adjuster (3) by turning counter-clockwise and slide the Adjuster directly below the Needle Handle (2) as shown in **Picture 2**. Turn the Needle Adjuster (3) clockwise to lock the Needle Handle (2)
- ! If any irregularities are noted, the instrument must be disposed of and replaced prior to use.

### **Handle illustration showing markings, unlocked and locked positions**

Markings on the Sliding Sheath Handle (8) and on the Needle Handle (2) indicate if the Adjusters are unlocked to facilitate movement and adjustments, or if they are locked in a safe and locked position.


<p><b><u>Unlocked</u></b> Adjuster is positioned off the symbol. The unlocking function is performed by turning the Adjuster counter-clockwise (needle tip towards distal).</p>	<p><b><u>Locked</u></b> Adjuster is positioned on the symbol. The locking function is performed by turning the Adjuster clockwise (needle tip towards distal).</p>
<p><b>Picture 3</b></p>	

### **Introducing the SonoTip® TopGain into the Endoscope**

- Slowly advance the SonoTip® TopGain Sheath (10) into the endoscope's accessory channel and Luer-Lock the handle (9) onto the endoscope's accessory channel Luer-fitting.
- Adjust the Sheath (10) length to the endoscope's working length by unlocking the blue Sheath Adjustor (7) and advancing the sliding Sheath Handle (8) until the tip of the Sheath (10) reaches the desired position.
- Secure the sheath length by locking the Sheath Adjustor (7), as shown in **Picture 3**.

### **Advancing the Needle**

- Unlock the Needle Adjuster (3) and slowly advance the Needle Handle (2) to extend the needle out of the Sheath (10) and towards the desired target lesion.
- **Round Stylet (11)**  
Loosen the round stylet from the proximal Luer-Lock of the Needle Handle (2) and retract the stylet of approx. 5 mm.
- Continue to advance the needle (11) slowly into the target lesion while maintaining constant visualization via the ultrasonic field. The needle projection length can be set at a maximum projection depth by locking the Needle Adjuster (3). Picture 4 shows the correct adjustment for a needle projection depth of 2 cm.

**Picture 4**

- ! Avoid abrupt movements while puncturing and projecting the needle in order not to harm proximate tissue.
- ! If resistance is felt at any time during use, do not exert excessive force as this may result in injury to the patient, endoscope and/or instrument.
- With the needle tip fully embedded into the target lesion, close the Stylet **(1)** again. Remove the Nitinol Stylet **(1)** completely from the Needle **(11)** and set it aside for reintroduction.
- ! Remove the Stylet carefully to avoid the splashing of possible contaminated tissue/body fluids to minimize risk of infection risk
- Attach the prepared aspiration syringe to the proximal end of the SonoTip® Easy Gain Pro and luer lock it securely into place. Release the stopcock on the syringe ¼ turn (aligning it parallel with the syringe). This will produce the negative suction pressure needed to aspirate the cell tissue into the needle's tip and lumen.
- Gently advance/retract the Needle Handle **(2)** into the biopsy site to insure adequate sampling.
- ! **Caution:** Avoid removing the needle **(11)** from the biopsy site during biopsy. This may compromise the integrity of the sample and diagnosis of the biopsy.
- With the needle **(11)** embedded in the biopsy site, close the stopcock on the syringe ¼ turn.
- Fully retract the Needle Handle **(2)** so that the Needle **(11)** is removed from the target lesion. Securely house the needle within the Sheath **(10)** by locking the Needle Adjuster **(3)** in the position as shown in picture 2.
- Open the Luer- Connector **(9)** from the endoscope and gently remove the SonoTip® TopGain from the endoscope.
- Direct the needle tip onto a slide and unlock the syringe plunger to eject the biopsy from the needle **(11)** onto the slide.
- ! The presence of a cytologist during extra-luminal endoscopic ultrasound-guided FNB may reduce the number of fine needle biopsy passes required from each lesion.
- If necessary, reinsert the nitinol stylet **(1)** into the needle **(11)** to help eject any remaining biopsy onto the slide. A new syringe with heparinized saline may also be used to flush the needle **(11)** of any excess tissue.
- Once the desired amount of biopsy has been collected, dispose of the entire needle system and syringe according to legal requirements and your facility's infection control protocols.

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**Medi-Globe GmbH, Medi-Globe-Str. 1-5, 83101 Achenmühle - Germany**
Medi-Globe GmbH